
CERNIERA F.M.
F.M. Hinge

CE
0 4 7 6

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***Unilateral solutions
applicable in 4 sectors***

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
The F.M. hinge is a cushioning extracoronar attachment for unilateral prostheses which is provided with an elastic return device to the static position. The system is made up of two separate and removable parts:

- the patrix, which is the hinge rotation axis, is connected to the fixed prosthesis;
- the matrix includes the male axis rotating seat, containing an elastic return device which can be removed from the distal end of the hinge.

The **patrix** is available in the following versions: the AISI 316L steel, the palladium alloy and the gold/pt alloy patrix, the castable resin patrix (which can be straight or inclined), and the HSL ceramic alloy patrix, which is used for the overmelting or welding process before or after the ceramization, with noble and non noble alloys. **Check the melting range of the alloys and solders to be used.** In castable resin and noble alloy versions, the transverse pin has been reinforced.

A new AISI 316L steel patrix with overmelting retention has been made with a reinforced transverse pin.

The AISI 316L steel, the palladium alloy and the gold/pt alloy **matrix** is provided with an hexagonal end for an higher retention. Thanks to the round grooves which are placed in the hexagonal part of the matrix, it can be used for welding and skeleton gluing techniques. The internal parts, which must be removed before welding in order to avoid the execution of heat treatments, are made up of a cylindrical pin with a ring stop in order to prevent the displacement, a spring, which can be provided with an higher strength upon demand, and a screw in order to fix the internal parts.

The F.M. hinge, rotating around the male transverse axis, determines free-end saddle sinking in a critical point ranging from 2 and 3 mm. It must be the only movement made by the hinge, as crosswise movements could damage the product. That's why it is forbidden to finish or polish the rotation axis surface or the device rotating body seat.

The device requires a correct relation between the fixed and missing elements to be inserted in the removable saddle. There isn't any fixed pattern, but the cases concerning a single abutment element are not taken into consideration (at least 2).

Periodic base reconstructions of the mobile part are highly recommended, so that the return device can go on performing its task.

Being a device for independent use, it doesn't need a parallelometer mandrel, but it is possible to use the New Ancorvis universal mandrel for the positioning of the patrix.

Especially for those who make use of this device for the first time, it is better to use the parallelometer plate as a reference, making sure that the model is not horizontally positioned, but slightly inclined on the vestibular surface. This inclination helps the patient to insert and remove the mobile part, for an higher stability during chewing process.

The metal casing which is included in the packaging must be used during the packing phase in order to prevent the resin from blocking the hinge, so that the attachment removability, due to wear, does not cause any crosswise movement of the hinge.

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
Instructions to be followed by the patient

- Keep a good oral hygiene
- Give instructions on how to fit the prosthesis
- In order to avoid excessive stresses on the attachments, it is advisable to reconstruct the saddle base if necessary.

SCHEDA TECNICA LEGHE PREZIOSE IMPIEGATE - TECHNICAL SHEET OF THE PRECIOUS ALLOYS WHICH ARE USED

		HSLCERAMICA HSL Ceramic	ORO Gold	PALLADIATA Palladium Alloy
Norma di riferimento - Reference regulation	UNI EN ISO	1562	8891	8891
Colore - Colour		Bianco - White	Giallo-Yellow	Bianco-White
Contenuto Metalli - Metal content	Au	60,0%	58,8%	10,0%
	Pt	24,9%	1,6%	
	Pd	15,0%	3,2%	20,9%
	Ir	0,1%	0,01%	
	Ag		22,4%	58,0%
	Cu		12,7%	9,5%
	Zn		1,3%	1,5%
	Ru			0,1%
Durezza Vickers - Hardness Vickers: HV 5/30	ricotto - annealed	150	150	160
	temperato - hardened	230	275	200
Intervallo di fusione - Melting point		1350-1460 °C	870-920°C	970-1020°C
Carico di rottura - Ultimate tensile strenght: N/mm ²	ricotto - annealed	580	530	650
	temperato - hardened	810	820	840
Carico di snervamento - Yield stress: N/mm ²	ricotto - annealed	450	450	260
	temperato - hardened	720	760	390
Densità - Density		18,1 g/cm ³	13,9 g/cm ³	11,2 g/cm ³
Allungamento - Elongation %	ricotto - annealed	18	20	25
	temperato - hardened	12	4	5

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The HSL ceramic alloy matrix can be used for the overmelting or welding process - Check the melting range of the alloys and solders to be used.

SCHEDA TECNICA ACCIAIO AISI 316 - Technical sheet of AISI 316 L stainless steel alloys

Norma di riferimento - Reference regulation	UNI EN ISO	10088-1
Analisi - Analysis: C 0,017 - Ni 10,02 - Cr 16,83 - Mn 1,43 - Si 0,43 - Mo 2,0 - P 0,033 - S 0,028 - N 0,06		
Durezza Vickers - Hardness Vickers: HV 5/30		200
Intervallo di fusione - Melting point		1371-1398 °C

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Antiresidue castable resin parts guarantee the absence of waste during the casting in the oven. The resin products which are present in the oral cavity are made up of a wear-resistant biocompatible material with an elastic memory according to the kind of attachments.

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