
LINEA MINI AN

Mini AN attachments

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From the analysis of the extracoronar attachment manufactured by New Ancorvis, it is clear that versatility is one of its main peculiarities, to planning the prosthesis working design perspective, meeting patient's needs, working on the resilience.

The assembling phases are not very different from the standard procedures related to the other extracoronar attachments, but, in order to achieve successful results, the attachments must be parallelly positioned using the parallelometer mandrel which is suitable for **MINI AN** attachments.


For a good operation, the attachments must be combined with a functional anatomic milling designed by the dental technician.

The **patrix**, which can have two inclinations (straight or a inclination of 30° or more) is available in the following versions: the castable resin and the HSL ceramic alloy version, which is used for the overmelting or welding process before or after the ceramization, with noble and non noble alloys. **Check the melting range of the alloys and solders to be used.** Each patrix can be integrated with four kinds of matrix which are interchangeable one with the other, but which have a different retention system.

Thanks to the side grooves, the AISI 316/L steel **matrix** can be used for welding and skeleton gluing techniques.


The operation must be carried out using the duplication matrix during the execution of the skeleton.

Before welding, remove the internal parts in order to avoid the execution of heat treatments.

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
MECHANICAL RETENTION SYSTEM

The system, which is made up of metal parts, is composed of a ring nut, a spring and a pin which can be provided in a larger size upon demand for wear prevention.

The internal parts can be easily replaced removing the ring nut from its housing. Thanks to the ring nut, it is also possible to increase or decrease pin projection using a special screwdriver.

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COMBINED RETENTION SYSTEM

It is made up of a spring and a pin (which is also available in a larger size) and a wear-resistant plastic sheath with two different strengths (red: lower retention – blue: higher retention) which guarantees a smoother friction during the insertion until the pin is released after having reached its end of stroke position.

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
RESILIENT RETENTION SYSTEM

The friction system is based on wear-resistant plastic sheaths instead of metal mechanical parts and it is possible to choose between two different strengths (red: lower retention - blue: higher retention), passing from a resilient retention system to a combined retention system replacing the sheath and inserting the metal mechanical parts

INCLINED REDUCED RESILIENT RETENTION SYSTEM

The resilient version has been integrated with a small-sized matrix for the execution of the prosthetic cases requiring the use of MINI A.N. attachments in the front areas in order to achieve better aesthetic results considering the small size of the attachment and for the solution of the cases requiring the use of an inclined matrix for free end saddle bone reabsorptions.


Fig. A
Picture A


Fig. B
Picture B


In case of the systems requiring the presence of sheaths, a white working sheath is already fitted in the attachment, so that dental technicians can carry out laboratory operations without damaging the final sheaths, to be replaced only at the end of the work.

The insertion will be made using a special insertion tool, pushing the sheath towards the two side retentions which are present in the matrix.

Check the direction of the sheath, making sure that the coloured border is directed towards the matrix (picture A).


Turn the insertion tool 180 degrees and insert the smooth side into the matrix for a few seconds for the adjustment of the sheath (picture B).

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Resiliente
ridotta
Reduced
resilient

Codice - Code	Matrice- Matrix	Patrice - Patrix
4300/d	Inox / Stainless Steel	HSL Ceramica dritta/ Straight HSL ceramic
4300/i	Inox / Stainless Steel	HSL Ceramica inclinata / Inclined HSL ceramic
4500/d	Inox / Stainless Steel	Calcinabile dritta / Straight casteable
4500/i	Inox / Stainless Steel	Calcinabile inclinata / Inclined casteable

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
La ceratura diagnostica del caso – *Diagnostic waxing of every single case*

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The pictures have been made in co-operation with TECNO.DENT dental laboratory, which gives theoretical and practical courses on the attachment assembly and milling. For further information, visit our website at www.newancorvis.com or call the toll free 800 850493 in order to receive the CD free of charge.

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Rigid anchoring systems are used for the construction of removable bridges and partial and combined prostheses (fixed and mobile prostheses) and, thanks to the technical and functional characteristics of the connection without resilience between the matrix and the patrix, the load is mainly distributed on the periodontium and the abutment teeth, guaranteeing prosthesis stability and functionality. Rigid anchoring systems are made up of primary crowns supporting the patrix and the milling, which is necessary to protect the attachment, the abutment teeth for the execution, the periodontium and consequently the device, and the removable part includes the matrix, the part which will fit into milling and the support for missing teeth, thus achieving a strong, functional and healthy prosthetic unit.

Instructions to be followed by the patient

- Keep a good oral hygiene
- Give instructions on how to fit the prosthesis
- In order to avoid excessive stresses on the attachments, it is advisable to reconstruct the saddle base if necessary.

SCHEDA TECNICA LEGHE PREZIOSE IMPIEGATE - TECHNICAL SHEET OF THE PRECIOUS ALLOYS WHICH ARE USED

		HSLCERAMICA HSL Ceramic	ORO Gold	PALLADIATA Palladium Alloy
Norma di riferimento - Reference regulation	UNI EN ISO	1562	8891	8891
Colore - Colour		Bianco - White	Giallo-Yellow	Bianco-White
Contenuto Metalli - Metal content				
	Au	60,0%	58,8%	10,0%
	Pt	24,9%	1,6%	
	Pd	15,0%	3,2%	20,9%
	Ir	0,1%	0,01%	
	Ag		22,4%	58,0%
	Cu		12,7%	9,5%
	Zn		1,3%	1,5%
	Ru			0,1%
Durezza Vickers - Hardness Vickers: HV 5/30	ricotto - annealed	150	150	160
	temperato - hardened	230	275	200
Intervallo di fusione - Melting point		1350-1460 °C	870-920°C	970-1020°C
Carico di rottura - Ultimate tensile strenght: N/mm²	ricotto - annealed	580	530	650
	temperato - hardened	810	820	840
Carico di snervamento - Yield stress: N/mm²	ricotto - annealed	450	450	260
	temperato - hardened	720	760	390
Densità - Density		18,1 g/cm³	13,9 g/cm³	11,2 g/cm³
Allungamento - Elongation %	ricotto - annealed	18	20	25
	temperato - hardened	12	4	5

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Antiresidue castable resin parts guarantee the absence of waste during the casting in the oven. The resin products which are present in the oral cavity are made up of a wear-resistant biocompatible material with an elastic memory according to the kind of attachments.