mega ***LIGHTMODEL***

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For further Informations contact us !

mega **LIGHTMODEL**

Light-curing and burn-out free of residue
modelling acrylic

mega **LIGHTMODEL**

Light-cured, ready-to-use material for modelling inlays, onlays, veneers, abutments and other modellable constructions as well as for blocking out undercuts.

Can be used in syringe (tup), ready for use

How to use:

mega LIGHTMODEL is ready to use and is taken directly from the dosing syringe. **mega LIGHTMODEL** can be adapted directly using the discharge cannula or applied and modelled with a modelling instrument.

mega LIGHTMODEL can be easily removed from the rotary syringe and adapted with the modelling instrument. After adaptation and shaping, **mega LIGHTMODEL** is light-cured in commercially available UVA and blue light polymerisation units (spectral range 320-500 nm). The curing time depends on the intensity of the light unit and can be 3-5 minutes. The modelled part is solid after light curing and can be processed immediately. The layer thickness should not exceed 3 mm per light curing. To fix soldered areas, the bridge parts are applied to the model and joined together at the separation points with **mega LIGHTMODEL**. After the light polymerisation of **mega LIGHTMODEL**, the joined bridge is lifted off the model. A soldering block is made with the soldering compound. By heating (preheating), **mega LIGHTMODEL** is expelled from the solder block. **mega LIGHTMODEL** burns without leaving any residue. The embedding technique is the same as for wax objects. Casting pins are heated and pressed or waxed into the hardened material. The casting muffle is burnt out after reaching a temperature of approx. 480°C, as **mega LIGHTMODEL** burns without leaving any residue.

Storage:

Store in a dark place and not above 25 °C / 77°F!

Precautions:

Do not expose the unprotected material to light before processing as it will polymerise. For professional use only, keep out of the reach of children.

durability:

24 month

Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

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