Heizplatte / Heating-Plate

Made in Germany

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Heating-Plate

Item-no. 100-466

Before start up please read the operating instructions of the equipment attentively.

Our **Heating-Plate** enables the quick and controlled pressure polymerisation of all auto-polymerisates.

The temperature of the heating plate is infinitely variable, so that the water for the polymerisation will be heated on max. 50° C, when using common polymerisation-pressure-pots. At the same time the control lamp shows the respective phase of heating. The 18 cm-diameter of the plate, allows the use of all common polymerisation pressure-pots.

Information

When using the Heating-Plate for the very first time, it is necessary to heat the appliance for approx. 10 minutes at maximum control-setting, without pot. This way, the surface of the heating plate gets the required protective effect.

A faster heating of the pot-content can be achieved by using warm water.

The Heating-Plate is equipped with a temperature-thermostat, so that it is impossible to overheat.

Guarantee

There is a 6 months guarantee on the Heating-Plate.

Safety data

A repair or an exchange of component parts should only be made by us or by traders, who are authorised by us. Also, the replacements of spare parts can only be carried out here in-house and only with original spare parts.

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