mega *PINKSIL N22*

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For further Informations contact us !

mega *PINKSIL N22*

addition-curing 1:1 duplicating silicone

mega-PINKSIL N22

is an addition-curing 1:1 duplicating silicone with final hardness of Shore A 22 and low viscous with best mechanical characteristics. It offers an extremely high accuracy in dimension and design of the duplicating form and the pink tint supports an optimal surface characterisation. **mega-PINKSIL N22** is very flexible in the final consistency despite high ultimate tensile strength and therefore easy to deflask particularly in case of undercuts. It is from extraordinary productivity since it contains no additional fillers and requires no wetting agents. Due to a platinum catalyst in component A the hardening without any shrinkage is ensured.

mega-PINKSIL N22 is also suitable for the casting without flasks and causes no limitation to the expansion value of the investment.

Working instructions:

mega-PINKSIL N22 component A and B must be mixed at room temperature in proportion 1:1 and filled into the dry plaster model containing flask constantly in the same place until the model is all round covered.

After a curing period of approx. 30 minutes **mega-PINKSIL N22** at a room temperature of 20-23°C the model can be replaced and should be left another few minutes to release probably existing reproductions of undercuts.

For further questions don't hesitate to contact us.

Technical data:

Components

	Component A	Component B
Viscosity:	2.500 mPa s	2.500 mPa s
Density at 20°C, DIN 51757:	1,1 g/cm ³	1,1 g/cm ³
Colour:	white	pink
Processing time:	5-6 min.	
Setting time:	15-30 min	
Colour:	pink	
Shore A after 30 min:	16-22, DIN 53 505	
Shore A after 60 min:	18-22, DIN 53 505	
Tear resistance:	3,5 N/mm ² , DIN 53 504 S 1	
Elongation at tear:	400 %, DIN 53 504 S 1	
Continuing tear resistance:	4,5 N/mm, ASTM D 624 B	
Linear change of measurement:	< 0,2 %	

These data represent approximate values and are not intended to produce specifications.

Storage durability:

mega-PINKSIL N22 is storable in original locked containers at least 12 months at a temperature between 10°C and 25°C.

Consider safety-related references of the safety data sheet

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