---

## maruvest *SPEED*

---

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# maruvest SPEED

Precision Shock Heat Investment  
for hot furnace

## maruvest-SPEED

is the latest development from **megadental GmbH**. This silica-phosphate bound, carbon-free investment can be used in both the shock heat technique and the conventional cold furnace technique. **maruvest-SPEED** is suitable for all Dental alloys. Due to the very fine grain structure, gives clean, bright castings which are precise and easy to de-vest. **maruvest-SPEED** offers high thermal stability and problem-free castings.

**Mixing ratios:** 13-14 ml liquid : 60 g **maruvest-SPEED** powder (sachet)  
36-38 ml liquid : 160 g **maruvest-SPEED** powder (sachet)  
22-24 ml liquid : 100 g **maruvest-SPEED** powder

Add powder to liquid and mix under vacuum for **75-90 seconds !!!** Allow to bench set for 30 minutes, Place in a cold furnace, burnout at 270-280°C and hold for 20-30 minutes (according to ring size), Hold the final end temperature for 30-60 minutes (according to ring size)

**Hot furnace technique:** After 12-15 min (**the first contact powder & liquid = 0 min. !!!**) put muffle in hot furnace (**by use without metal-ring – 820°C**) or (**by use with metal-ring max. 850°C**), hold at this temp. for 15 min. before going up to the final temperature and holding for 45-60 min. according ring size before casting.

**Cold furnace** Place the investment ring into either a cold furnace or one that is at max. 150-200°C and raise **technique:** to 270-280°C to burnout wax for 20 minutes, than up to 570-580°C and hold for 20-30 minutes after which increase to final end temperature and hold for 30-45 min. according ring size before casting.

**De-Vesting:** Allow the ring to cool slowly, do not quench as this can result in the oxide layer adhering to the alloy.  
**Casting temperatures:** Precious and Semi-Precious alloys: 700°C - 850°C  
Non-Precious: 860°C - 950°C  
Please also observe the instructions of the alloy manufacturer

**Storage:** can be stored in a dry and frost-free room for up to 8 months,

**Caution:** Investments contain quartz-wear adequate protection-do not breath investment dust!!!  
If the maruvest Liquid comes into contact with skin wash off with soap and water. If **maruvest liquid** comes into contact with eyes wash out immediately and seek medical attention !!! **maruvest Liquid** is alkaline. Consult **maruvest-SPEED** safety data sheet for further information

### Expansion properties:

Liquid-%	Setting-Expansion	Thermal-Exp. (750°C)	Total-Exp. (750°C)
Red Special (100%)	2,0-2,2%	1,2 %	3,2-3,4 %
Blue Regular (100 %)	1,6-1,8 %	1,2 %	2,8-3,0 %
Red Special (75%)	1,2-1,6 %	1,2%	2,4-2,8 %
Blue Regular (75 %)	0,7-1,0 %	1,2 %	1,9-2,2 %
Red Special (50%)	0,7-0,8%	1,1 %	1,8-1,9 %
Blue Regular (50 %)	0,5-0,6 %	1,1 %	1,6-1,7 %
Graphit Liq. (100 %)	1,0-1,2 %	1,2 %	2,2-2,4 %

### Suggested Liquid Concentrations:

		Blue Regular Liq.	Red Special Liq. GRAPHITE Liq
<b>Precious Alloys</b> Inlay, Crowns, Bridges...	65 - 85 %	60 - 70%	80 - 100 % as Gold part...
<b>Non-Precious alloys</b> Crowns, Bridges, Telescopic works...	95 - 100 %	90 - 100 %	Not usable ! Do not use with Palladium

All figures are approximate. When diluting **Liquid**, than use only distilled water.  
The above are given as a guide only and may be varied according to individual technique.

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