maruvest *Cerampress*

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Precision Shock Heat Investment for hot furnace
indication for pressing ceramic

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maruvest Cerampress

Precision Shock Heat Investment (650°C-950°C) for hot furnace indication for pressing ceramic

maruvest-Cerampress

is the latest development from **megadental GmbH**. This silica-phosphate bound, carbon-free investment can be used in both the shock heat technique and the conventional cold furnace technique. **maruvest-Cerampress** have the indication to use by pressing ceramic. Due to the very fine grain structure, gives clean, bright pressing which are precise and easy to devest. **maruvest-Cerampress** offers high thermal stability and problem-free use.

Mixing ratios: 23-24 ml liquid : 100 g **maruvest-Cerampress** powder (sachet)

Add powder to liquid and mix under vacuum for **75-90 seconds** Allow to bench set for 60 minutes. Place in a cold furnace, burnout at 270-280°C and hold for 20-30 minutes (according to ring size), Hold the final end temperature for 30-60 minutes (according to ring size). Hot furnace. After 12-15 min (the first contact powder & liquid = 0 min.) put muffle in hot furnace technique: **((by use without metal-ring – 800°C) or (by use with metal-ring max. 850°C))**, hold at this temp. for 15 min. before going up to the final temperature and holding for 45-60 min. according ring size before pressing.

Cold furnace technique: Place the investment ring into either a cold furnace or one that is at max. 150-200°C and raise to 270-280°C to burnout wax for 20 minutes, than up to 570-580°C and hold for 20-30 minutes after which increase to final end temperature and hold for 30-45 min. according ring size before pressing.

De-Vesting: Allow the ring to cool slowly, do not quench !

Suggested Use instruction from ceramic manufacturer.

Casting temperatures: Please note the details of the ceramic manufacturer.

Storage: can be stored in a dry and frost-free room for up to 8 months,

Caution: Investments contain quartz-wear adequate protection - do not breath investment dust!!!

If the **maruvest** Liquid comes into contact with skin wash off with soap and water. If **maruvest** liquid comes into contact with eyes wash out immediately and seek medical attention !!! **maruvest** Liquid is alkaline. Consult **maruvest-Cerampress** safety data sheet for further information.

Expansion properties:

Cerampress-Liquid-%	Setting-Expansion	Thermal-Exp. (750°C)	Total-Exp. (750°C)
(100 %)	1,2 %	1,2 %	2,4 %
(80%)	0,7 %	1,2%	1,9 %

All figures are approximate. When diluting **Cerampress Liquid**, than use only distilled water. The above are given as a guide only and may be varied according to individual technique.

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