mega *CREATIV*

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For further Informations contact us !

mega *CREATIV*

High pigmented acrylic colors
Autopolymerizable, MMA-based Powder

mega *CREATIV* allows individual customization of the artificial gingiva of dentures. It can be used while preparing new dentures or for additional coloring of already completed dentures. It is perfectly suitable for the **mega CAST** System in connection with **megaCRYL N**, our unique micro pearl denture base resin.

mega *CREATIV* is based on MMA. The intensive colors are color stable, free from tertiary amine and contain organic and inorganic, biologically compatible color pigments.

Use: The application of the colors should be carried out after all other necessary preparations such as boiling-out of wax, conditioning and repositioning of the teeth and isolation of the model with **megasol NF** have been fulfilled and the flask has been prepared for the casting process. **megaCRYL** Monomers have to be used. The colors are miscible with each other and can be intermixed with anyone of our denture base resins. **mega *CREATIV* Colors** are applied with the brush-on technique by moistening the tip of the brush with monomer and then taking up a bit of powder until a thixotrope drop is left on the tip of the brush which is then placed inside the flask. The material should be stable before being applied inside the flask to ensure that it does not run into areas where it is not wanted. Then, the material can be applied drop by drop. (use appropriate, monomer resistant brush).

Important: Once the material has started to dry (may result in a whitish coloration) immediately moisten the surface with monomer. It is suggested to cover the entire preparation after application of the colors with monomer before closing the flask and casting the denture base resin. Polymerization is carried out as usual, preferably in a pressure polymerization device.

mega *CREATIV* can be combined with any denture base resin.

Homogeneity:

Combining mega *CREATIV* Colors with megaCRYL S+N results in a one-piece denture, consistent of a single material in different individual shades. The denture does not show streaks or joints between the colors and can be grinded and polished as usual.

Delivery form: mega *CREATIV* Set with 8 colors containing 8g of powder. Each color is also sold separately.

Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

megadental is certified in accordance to ISO 9001 and guarantees efficient quality through maintaining a deliberate quality management system. Product properties as mentioned in this instruction can only be guaranteed when the procedural steps are carried out properly. The user is responsible for correct handling of the product. Any defective results cannot be accounted for, because megadental has no influence on the subsequent processing by the user. Any claims for damages can solely be granted for the value of our goods.

according to medical device directive 93/42 EEC annex IX class **Ila** for fixed and removable dentures
according to DIN EN ISO 20795-1 Type 2 Group 1

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