megaCRYL S + N

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megaCRYL S + N

Auto polymerizing denture base acrylics
with micropearls

megaCRYL S + N

are cold curing acrylics with micropearls based on methylmethacrylate. The materials are suitable for any indication. They are free of cadmium and tertiary amines and therefore **color stable**. **megaCRYL S + N** are easy to use and they are characterized by their high mechanical values, natural shading as well as fast and safe polymerization.

Preparation: The arrangement of the prosthetic teeth in wax need to be fixed with a key made of plaster or silicone. Key, model and teeth are boiled out thoroughly to remove any wax residue. The plaster model needs to be watered for at least 10 minutes to assure that no air is pressed from the plaster into the acrylic or the isolation. To provide isolation between the acrylic and the plaster alginate based isolation is used. **Important ! Two thin layers of the isolation need to be applied ! Make sure the first layer of isolation is not completely dried before applying the second layer for this may cause white marks on the acrylic ! Before the isolation gets in contact with the acrylic both layers need to be completely dried !**

Mixing: **Suggested mixing ratio: 10 g of powder with 4 - 5 g of liquid.**
30 g powder with 13 - 15 g of liquid. (for a normal size OK-denture)

Place required amount of liquid into the mixing cup (for example **megaEXACT-kit**), then add the according amount of powder. **Important !** Because of the micropearls it is essential to assure a "lump-free" and homogenous mixture.

After a short thickening phase of approx 20-30 seconds the material should be stirred one more time ! Dosing individually is possible, yet weighing of the components is highly recommended ! A higher amount of liquid results in a greater shrinkage !

Processing times: **Pouring phase:** **megaCRYL N** 2-3 minutes, material then becomes stable.
megaCRYL S 1-2 minutes, material then becomes stable.
Plastic modeling phase: **megaCRYL N** 3-5 minutes.
megaCRYL S 2-4 minutes.

The provided processing times refer to a material- and room temperature of 23°C [73,5°F]. Higher temperatures extend, lower temperatures shorten these processing times.

Processing: The mixed material is poured into the precasts / injection cylinders within the pouring phase. As soon as the plastic-modeling phase begins the material becomes stable, stops flowing out of the precasts and can be modeled. The **polymerization** is to be carried out after a maximum of **6-8 minutes** in a pressure polymerization unit for 20-30 minutes at a pressure of 2-3 bar and a water temperature of 40-50°C.

When used in the flask casting technique, e.g. the mega-CAST System and using agar-based duplicating gel or mega-Pinksil 22 duplicating silicone, the polymerization time is 30-40 minutes.

Colours: pink, pink-veined, light-pink, light-pink-veined, pink-opak, V5, clear

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Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

Warranty: All references are based on our present knowledge. The user himself is responsible for the correct use and safe handling. Should any claims for indemnification be asserted, they shall only be applicable for the value of the delivered goods.

according to medical device directive 93/42 EEC annex IX **class IIa** for fixed and removable dentures
according to DIN EN ISO 20795-1 Type 2 Group 1

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