mega PUR

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mega PUR

Polyurethane-based, two-component die material for high-precision models.

mega PUR is suitable for fabricating precision-, master, saw- and special models. Especially when making fragile preparations or very thin layered posts, very high details and dimensionally stable results are achieved. The material is easy to handle. It is designed to make removable gum preparations in implantology and suitable for both partial or complete models.

mega PUR is designed for use with silicone impression materials. Before using polyether or polysulphide impression materials (e.g. Impregum or Permedyne) mega PUR separation spray should be applied to the impression. Apply the separation spray by thinly spraying from a distance of approx. 20-30 cm and avoid forming pools. Allow 3 minutes of air drying. **mega PUR** is not suitable for hydrocolloids and alginate impression materials.

mega PUR is fluid and easy to pour (even with very intricate cavities), produces precise details and is dimensionally stable. It also guarantees optimum edge stability and abrasion resistance. It is resistant to pressure, impact and shear forces and can easily be trimmed and sawed. When used correctly, no loss of structure due to washing, steam cleaning or ultrasonic cleaning can be guaranteed.

Instructions for use:

1. Thoroughly stir the base material with the wooden spatula before use.
2. Remove the desired amount of base material and place into the measuring beaker.
3. Add the according amount of hardener.
4. Stir thoroughly for at least 25 seconds until both components become a homogenous mass with no visible streaks.
5. Pour the die material thinly into the cleaned (separated) and dried impression or silicone cast. The working time is approx. 3-4 minutes at 20-23°C.
6. To prevent air bubbles, the impression can be placed on a vibrato at low setting when pouring and/or cured for approx. 10-15 minutes at approx. 2-4 bar in a dry pressure chamber.
7. The material is set after 60 minutes and the impression can be removed.
8. The model can be trimmed, ground and drilled after approx. 120-180 minutes.
9. mega PUR is compatible to all common die spacer lacquers and is resistant to separating agents.

Mixing ratio:

(parts by weight: 100 component A : 20 component B)

Base / g	Hardener / g
10 g	2 g
15 g	3 g
20 g	4g
30 g	6 g
40 g	8 g

Color:

Component A:	white
Component B:	brown
Working time:	3-4 minutes
Removable after:	60 minutes
Setting time:	2-3 hours
Shore-D (24h):	72
Bending strength (7dRT/24h 100°C) :	approx. 65 N/qmm
Contraction:	approx. 0,03 %

For a medium row of teeth you need approx. 20-25 ml of material

Important notice:

The material is extremely sensitive to moisture. Avoid contact of the components with moist or wet surfaces. Keep containers closed immediately after use. Isopropanol or alcohol is suitable for cleaning the working utensils. Do not expose to heat, the material sets faster at higher temperatures. Additional information is provided in the Material Safety Datasheets.

Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

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