mega RESIN X-pose

CT & DVT

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For further Informations contact us !

mega *RESIN X-pose*

cold curing, radio-opaque acrylic

mega-RESIN X – pose CT

is a self curing acrylic for the manufacture of removable diagnostic aids, planning teeth and scan dentures. This acrylic is perfectly suitable for 3D X-ray diagnostics in CT and the following prosthetic backward scheduling. The perfectly balanced radiopacity guarantees clean determination of the prosthetics in all standard software applications.

mega RESIN X – pose DVT is available in two different radiolucency concentrations. In low-radiation recordings with DVT machines, mega Resin X-pose DVT shows no so-called artefacts. The conversion and segmentation of DVT – DICOM data is simplified and can therefore be carried out even more economic.

mega RESIN X - pose shows excellent grinding and polishing features after polymerization.

Preparation: The plaster model should be watered for at least 10 minutes so that no air is discharged from the plaster into the acrylic or the isolation which would then result in bubbles. Always use alginate based isolation fluids. (Plaster against acrylic) **Important!** Apply two thin layers of isolation fluid. Make sure the first layer is not completely dried before applying the second layer. Before getting in contact with the acrylic, the second isolation layer should be completely dry.

Mixing: Mixing ratio (by weight) **10 g of Powder and 4-5 g of liquid**. Dosing individually is possible, when specific flow-abilities are needed. When the material has reached the desired consistency, the material can be handled as usual.

Processing: After the swelling phase the acrylic is filled into the precast. As soon as the materials plastic phase begins, the material stops running and will not flow from the precast and can be modeled.

The polymerization is carried out in a pressure polymerization unit in warm water (37-45 °C) and a pressure of 2-4 bar. The polymerization time is 20-30 minutes. Finishing and polishing is carried out as usual.

Dental technician's tip:

pre-warming the cast (40°-50°C) and exact dosing of the components assures a perfect polymerization result.

Suggestion: It is recommended to store the device in warm water (35-38°C) for at least 24h in order to prevent irritations of taste due to residual monomer.

Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

Classification according to MDD 93/42/EEC Annex IX **Class I** for auxiliary dental materials with short-term intraoral use

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