mega *LIGHTTRAY*

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For further Informations contact us !

mega *LIGHTTRAY*

Light curing resin for fabrication of custom trays
pink, blue and transparent

How to use:

The edge course for the individual spoon is marked on the plaster model. Attach space-holders if required. If no place-holders are needed, insulated the plaster model with **mega *LIGHTTRAY***.

Remove protective film and place the **mega *LIGHTTRAY*** on the model, shape it and cut at the edges as marked. The transparency of the **mega *LIGHTTRAY*** transparent allows very accurate shaping. Excess material can easily be used to make the handle.

is polymerised in a normal light device (blue or UVA light). Illuminate the shaped spoon for about 5 minutes from behind.

The shaped spoon is ready for use and further processing once it has been polymerised. Any slight inhibition layer produced by oxygen in the air can be removed by wiping the spoon with a cloth dipped in **mega *Multiclean*** or an alcohol-based detergent.

Storage:

Store in a dark place and not above 25 °C / 68°F!

Precautions:

In its unhardened state contains methacrylic acid ester, which can cause skin allergies. Operators who are allergic or hypersensitive should use appropriate protection. It is essential that polymerisation is sufficient and complete, as illumination for an insufficient length of time or at insufficient intensity leads to a loss of strength.

durability:

24 month

Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

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