---

## **megafeel-EXACT**

---

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***For further Informations contact us !***

---

# megafeel-EXACT

## Reversible Duplicating Gel

**megafeel-EXACT** is a reversible duplicating gel, especially designed for the casting technique. It's unique feature is the low working temperature of 40°C – 43°C. Therefore the difference between the working temperature and the room temperature is just marginal. This guarantees a very stable and precise consistency of the duplicating mold.

### Processing:

**Preparation:** The Gel should always be cut into small pieces before being placed into the melting device. It should then be heated whilst **stirring constantly** up to a **maximum temperature of 90°C-95°C**, until reaching a homogenous consistency. The duplicating gel must not boil, for that drastically reduces the materials product life. **In order to prevent the material from dehydrating always keep the lid of the melting device closed.** The gel should then be allowed to cool down to a temperature of **40°C - 43°C**. Keep stirring constantly while cooling down and hold the appropriate working temperature.

**Information:** **When melting the gel for the first time no water should be added.** After use, the duplicating models should always be „washed“ with water and stored in a sealed container until being re-used. The „washing“ thereby compensates the loss of water. But still it is recommended to add about 50 ml - 100 ml of water (per 6kg of gel) after every third melting procedure. **Avoid adding too much water because that would lower the materials viscosity.** Also avoid excessive contact of the gel with acrylic monomer!

**Moisten the models prior to the duplication process for at least 10 - 15 minutes in hand warm water!**

### Product specifications:

**Melting temperature:** 90°C - 95°C

**Working temperature:** 40°C - 43°C

**Re-Usability:** 5-10 times (depending on working conditions)

**Shelf life:** 4 years when stored in original container between 5°C and 26°C

### Package:

**megafeel-EXACT:** 6 kg PE-Container, Turquoise or nature  
Please always store used and washed gel in closed original container!

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