---

## **megafeel**

---

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***For further Informations contact us !***

---

# megafeel

Reversible duplication gel

**megafeel** a reversible duplicating gel on base of agar – agar in green, blue and pink colour suitable for all gypsum-free investment materials (**megafeel**), for mega-CAST acrylic casting technique (**megafeel Extrahart**), and the special duplicating gel (**megafeel GIPS**) for plaster models.

Due to its excellent reproduction accuracy even to the smallest details **megafeel** is very well suitable for demanding partial denture. It produces a tear proof, highly elastic duplicating mould which doesn't cause any problems either at undercuts.

A high-quality base material produces an extremely long durability.

**Working instructions:** Reduce the duplicating gel to pieces as small as possible and fill it into the duplicating unit which is to heat up to **90°C – 95°C** under constant stirring.

**Don't add water with first melting.** Constant stirring is important to avoid lump formation and to get the gel well-balanced liquid.

An overheating up to cooking should be absolutely avoided as this can belittle the durability. In order to exclude unnecessary evaporation from liquid, keep the equipment always closed.

After melting the gel must be cooled down to:

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This working temperature should be kept under constant further stirring.

Before duplicating the model must be put **15 – 20 min. into slightly warm water.**

The duplicate must be absolutely dry when pouring in the gel to avoid the sticking of investment at the duplicate. Bonding time after duplicating at room temperature is approx. 1 h.

Wash the duplicate thoroughly with water after took out of the investment model and keep it in a closed container. The loss of water is replaced as far as possible by washing. However before melting again an addition of 50–100 ml can be necessary if the gel became too dry during the storage.

This addition of water should be limited to the absolute necessary quantity, in order not to thin down the duplicating gel excessively.

When using the duplicating gel in acrylic prosthetics it's important to permit contact with monomer only in inevitable degree.

#### Packing:

<b>megafeel:</b>	6 kg PE- bucket, green.
<b>megafeel Extrahart:</b>	6 kg PE- bucket, blue and transparent.
<b>megafeel Gips:</b>	6 kg PE- bucket, pink.

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