megamix *SPEED*

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For further Informations contact us !

megamix SPEED

Precision Shock Heat Investment
for hot furnace.

megamix SPEED

is the latest development from **megadental GmbH**. This silica-phosphate bound, carbon-free investment can be used in both the shock heat technique and the conventional cold furnace technique. **megamix SPEED** is suitable for all Dental alloys. Due to the very fine grain structure, gives clean, bright castings which are precise and easy to devest. **megamix SPEED** offers high thermal stability and problem-free castings.

Mixing ratios: 22-23 ml megamix SPEED Liquid : 100 g megamix SPEED Powder
35 ml (+/- 1 ml) megamix SPEED Liquid : 150 g megamix SPEED Powder
106 ml (+/- 2 ml) megamix SPEED Liquid : 450 g megamix SPEED Powder

Add powder to liquid and mix under vacuum for **75-90 seconds !!!** Allow to bench set for 30 minutes, Place in a cold furnace, burnout at 270-280°C
Hold for 20-30 minutes (according to ring size), hold the final end temperature for 30-60 minutes (according to ring size)

Hot furnace technique: After 12-15 min (**after mixing powder & liquid**) put muffle in hot furnace (**by use without metal-ring max. 820°C**) or (**by use with metal-ring max. 850°C**), hold at this temp. for 15 min. before going up to the final temperature and holding for 45-60 min. according ring size before casting.

Cold furnace technique: Place the investment ring into either a cold furnace or one that is at max. 150-200°C and raise to 270-280°C to burnout wax for 20 minutes, than heat slowly up (5-7°C/min.) to 570-580°C and hold for 20-30 minutes after which increase to final end temperature and hold for 30-45 min. according ring size before casting.

De-Vesting: Casting Allow the ring to cool slowly, do not quench as this can result in the oxide layer adhering to the alloy.
Precious and Semi-Precious alloys: 700°C - 850°C

temperatures: Non-Precious: 860°C - 950°C
Please also observe the instructions of the alloy manufacturer

Storage: can be stored in a dry and frost-free room.

Caution: Investments contain quartz. Do not breath investment dust!!! Please use appropriate respiratory protection. If the **megamix SPEED** Liquid comes into contact with skin wash off with soap and water. If **megamix liquid** comes into contact with eyes wash out immediately and seek medical attention !!! **megamix** Liquid is alkaline. Consult **megamix SPEED** safety data sheet for further information

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Suggested Liquid Concentrations:

for the model	55 % - 75 %
for the muffle / cylinder	30 % - 40 %

All figures are approximate. When diluting **Liquid**, than use only distilled water.
The above are given as a guide only and may be varied according to individual technique.

Important Note: Do not bring **megamix SPEED** with gypsum-based materials in connection.
Storage stability: **megamix SPEED** stored in dry and frost-free areas

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