megamix

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megamix

Model Casting Investment
for all Co.-Cr.-Alloys

megamix

is a precision model casting investment for all Co.-Cr.-alloys, high accuracy, simple and easy handling, dependable in use. **megamix** can be used with **mega-Cristall** expansion liquid or with water. **megamix** is specially tested with the duplicating **megafeel** and the cold dipping hardener **mercapor**.

Mixing ratio: 14 ml-15 ml liquid / 100 g **megamix** powder
56 ml-60 ml liquid / 400 g **megamix** powder

Working instructions: Please use **megamix** by temperatures 18°C-22°C The working time is 4-5 minutes. Add powder to liquid and mix under vacuum for 60-75 sec.

Model: Setting time for the model 25-30 minutes by **megafeel**, by silicon 15-20 minutes.

Drying and Hardening the model: 45-60 minutes at a temperature of between 200°C-220°C in the dryer. models ca. 2-3 seconds in a cold dipping hardener mercapor , if you use silicon duplicating you don't need the drying and hardening.

Muffles: Use **Water** for the master. The use of metal rings is not required! If you work with metal rings, please use generally 2 fleece deposits. The fleece deposits should not be dry (use Vaseline) and please leave overlap.

Setting and preheating time: After 15-20 minutes setting-time, place in a cold furnace or (max. 200 C°) warm furnace burnout at 270-285°C, hold for 20-30 minutes, then burnout at 400°C and hold, hold for 20-30 minutes then burnout at 580°C and hold for 20-30 minutes. then go to final end temperature, hold for 45-60 minutes (according to ring size).

Final end-temperature: Partial plates 880°C-930°C, total plates 920°C-950°C . Please also observe the instructions of the alloy manufacturer.

Storage: **megamix** can be stored in a dry and frost free room for up to 8 months.

Caution: *Investments contain quartz, wear adequate protection do not breath in dust. megamix should never get into contact with plaster or plaster containing materials.*

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