---

# meganium CS

---

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***For further Informations contact us !***

# meganium CS

Nickel-chrome metal-to-ceramic alloy

## Chemical analysis:

63 % Ni  
23 % Cr  
10 % Mo  
1,5 % Fe  
2,0 % Si  
0,05 % C  
<1 % Ce

## Alloy characteristics:

hardness: HV 190 – 200  
elongation: 12-15 %  
Tensile Strength: 450-550 N/mm<sup>2</sup>  
density: 8,2 g/cm<sup>3</sup>  
expansion-cof.: 14,0 x 10<sup>-6</sup> (20-600°C)  
melting range: 1250 – 1310°C  
casting-temp.: 1410 – 1440°C

(Specification of average value according to reference-analysis)

CE 0297 / class IIa / Medical Devices Directorate 93/42/EEC

**Material:** meganium CS is a Ni.Cr.Mo- metal-to-ceramic alloy, containing no precious metals, whose excellent, finely-balanced properties considerably fulfils the values conform standard **DIN EN ISO 22674:2006, DIN EN ISO 9693.**

The use of materials largely free from impurities and special melting and manufacturing processes impart extreme purity to meganium CS. It is free from beryllium. Very easy and sure to work.

**Preheating:** Crown and bridges temperature 800 °C - 850 °C

**Melting/casting: Precept:** Do not overheat alloy. Use only clean crucibles. Have one crucible per alloy.

**Recommendation:** To make an exact identification of each batch possible, cast new metal, invariably.

**If re-casting:** Portion of scrap metal max. 50% and only of the same batch as the new metal. Never mix different batches!

Use ceramic crucibles without exception. Vacuum pressure casting with induction heating and centrifugal casting with induction heating.

When the last solid component has submerged completely in the melt, continue to heat for 0 to 12 seconds depending on the induction capacity of the casting unit, then trigger casting.

Flame centrifugal casting: Cast when the last solid component has submerged completely in the melt and the melt clearly moves due to the flame pressure.

**Ceramic:** Use ceramics according to DIN EN ISO 9693 with firing temperatures of up to approx. 980 °C (e. g. Carat, Biodent, Duceram, IPS-Classic, Omega, VMK 95).

Also suitable: Ceramics with reduced firing temperature (e. g. Omega 900). Always follow the ceramic manufacturer's instructions!

Always blast the surface to be veneered (pure aluminium-oxide 110 – 250 my 3 - 4 bar) and clean the frame thoroughly (steam clean or boil in aqua dest.). After cleaning, hold with artery clamps and refrain from touching.

**Oxide firing:** Not necessary. If oxide firing is carried out to check the surface (900 °C / 5 minutes), always blast the oxide again (with pure aluminium-oxide 110-250 my 3 - 4 bar).

**Final work:** Blast visible metal surfaces with pure aluminium-oxide, then blast-polish external surfaces with mega-Perlblast 50 – 125 my. Clean thoroughly (steam clean or boil in aqua dest.).

**safty hint:** Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask type FFP3-EN149:2001!

**Secondary effects:** Such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

**Reciprocal actions:** In case of occlusal or approximate contact of different alloys electrochemically based reactions may very rarely occur.

**Reactions:** In case of known incompatibilities and allergies to contents of the alloy.

**Warranty:** All references are based on our present knowledge. The user himself is responsible for the correct use and safe handling. Should any claims for indemnification be asserted, they shall only be applicable for the value of the delivered goods.

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**For further Informations contact us !**

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