meganium **NF**

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CE 0297 / class IIa / Medical Devices Directorate 93/42/EEC

Material: **meganium NF** is the Co-Cr-Mo dental alloy, without nickel and beryllium, whose excellent, finely-balanced properties considerably exceed the values conform standard **EN ISO 22674:2006**. The use of materials largely free from impurities and special melting and manufacturing processes impart extreme purity to **meganium NF**. It contains neither nickel nor beryllium as alloying elements and displays excellent tissue compatibility & very easy to work with.

Preheating: Crown and bridges temperature 950 °C - 1050 °C

Melting/casting: **Precept:** Do not overheat alloy. Use clean crucibles, only. Have one crucible per alloy.

Recommendation: To make an exact identification of each batch possible, cast new metal, invariably.
If re-casting: Re-cast identic alloys, only. Blast old material and add at least 50 % of new material. Use ceramic crucibles without exception. Vacuum pressure casting with induction heating and centrifugal casting with induction heating. When the last solid component has submerged completely in the melt, continue to heat for 0 to 12 seconds depending on the induction capacity of the casting unit, then trigger casting. Flame centrifugal casting: Cast when the last solid component has submerged completely in the melt and the melt clearly moves due to the flame pressure.

Final work: Blast visible metal surfaces with pure aluminium-oxide, then blast-polish external surfaces with mega-Perlblast 50 – 125 my. Clean thoroughly (steam clean or boil in aqua dest.).

safty hint: Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask type FFP3-EN149:2001!

Secondary effects: Such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

Reciprocal actions: In case of occlusal or approximate contact of different alloys electrochemically based reactions may very rarely occur.

Reactions: In case of known incompatibilities and allergies to contents of the alloy.

Warranty: All references are based on our present knowledge. The user himself is responsible for the correct use and safe handling. Should any claims for indemnification be asserted, they shall only be applicable for the value of the delivered goods.

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