mega Model Resin NF

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Rapid hardening model resin

mega-Model Resin

- for conical caps and variations of the telescope caps for secondary parts (conical caps, telescope caps, RSS-bed, latch constructions, anchor strap caps, etc.)
- for inlays, subligule clip and palatinal constructions, connection parts, clips
- for adhesive bridges
- as soldering fixing for the embedding

physical values

setting time (23°C),	4 minutes
bending resistance (37°C, 10 min.),	>650 Kgf/cm ²
setting contraction (23°C) 30 min.,	~ 0,5-0,65 %

instruments and material:

resin mixing cup, brush, tissues, etc.

powder and liquid preparation:

The amount of powder necessary is filled into one cup and drops of liquid into the other.

cap construction:

(a outside cap as an example – secondary part – on the conical inside cap – primary part -):

Moisten the top of the small brush with liquid and take up some powder with it. This mixture is applied on the inside cap. Repeat the process of material taking and removing until the cap is complete. An isolation material for the inside cap is not needed. Clean the brush with liquid and a tissue after finishing the work or when necessary during application. Do not leave any resin remains on the brush.

controlling the interior area of the cap:

When the resin is hardened by polymerization, take off the cap of the primary cap. Taking up by pressing the resin edge. The resin gets a mat look when removing. Now you can remove the material from the primary part undamaged.

precision work:

Put the cap again on the inside cap and grind it back to an even thickness of 0,3 to 0,4 mm with an appropriate carbide milling cutter. The precision work must be done with the lowest number of revs possible, because the resin is warmed up to much at high numbers of revs and therefore would cause distortions. Check the thickness of the cap with sensor compasses or watch.

wax modelling and facing window:

The final cap form is modelled with inlay wax. Then the window for facing and retention barriers are applied. By that only the surface could get rough.

embedding, pre-heating and founding:

These processes are carried out in the usual way. For a single cap (weight under 0,3 g) the pre-heating is done similar to the regular wax modelling.

For large bridges about one hour is needed for burning out at about 270°C before the founding temperature is set.

Do not use any tension reducing media.

Attention: Not suitable for use intraoral!

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For further Informations contact us !

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