

according to Regulation (EC) No 1907/2006

mega TRAY NF

Revision date: 14.05.2018

Product code: 21371_EN

UFI: W9RW-54JJ-P60G-NNJS

Page 1 of 7

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

mega TRAY NF Product group:

Powder component

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Dental appliance.

1.3. Details of the supplier of the safety data sheet

| Company name: | megadental GmbH |
|-----------------|----------------------|
| Street: | Seeweg 20 |
| Place: | D-63654 Büdingen |
| Telephone: | ++49 (0) 6042-9755-0 |
| e-mail: | info@megadental.de |
| Contact person: | Murat Büyük |
| e-mail: | info@megadental.de |
| Internet: | www.megadental.de |

Telefax: ++49 (0) 6042-9755-20

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Regulation (EC) No. 1272/2008

Hazard categories: Respiratory or skin sensitisation: Skin Sens. 1 Hazard Statements: May cause an allergic skin reaction.

2.2. Label elements

Regulation (EC) No. 1272/2008

Hazard components for labelling dibenzoyl peroxide; benzoyl peroxide

Signal word:

Pictograms:



Warning

Hazard statements

May cause an allergic skin reaction.

Precautionary statements

Avoid breathing dust/fume/gas/mist/vapours/spray. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

2.3. Other hazards

No data available

SECTION 3: Composition/information on ingredients

3.2. Mixtures



according to Regulation (EC) No 1907/2006

mega TRAY NF

Revision date: 14.05.2018

Product code: 21371_EN

Page 2 of 7

Chemical characterization

Mixture on base of acrylic powder and catalyst.

Hazardous components

| CAS No | Chemical name | | Quantity | |
|--|--------------------------------|-----------------------------|-----------|---------------|
| | EC No | Index No | REACH No | |
| | Classification according to F | Regulation (EC) No. 1272/20 | 008 [CLP] | |
| | Acrylic polymer | | | 95 - <= 100 % |
| | | | | |
| | | | | |
| 94-36-0 dibenzoyl peroxide; benzoyl peroxide | | 1 - < 5 % | | |
| | 202-327-6 | | | |
| | Org. Perox. B, Eye Irrit. 2, S | Skin Sens. 1; H241 H319 H3 | 17 | |

Full text of H and EUH statements: see section 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

Remove contaminated, saturated clothing immediately. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

After inhalation

Provide fresh air. If breathing is irregular or stopped, administer artificial respiration. Medical treatment necessary.

After contact with skin

After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

Rinse mouth immediately and drink plenty of water.

4.2. Most important symptoms and effects, both acute and delayed

No information available.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

Unsuitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixture

Non-flammable.

5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

Additional information

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

according to Regulation (EC) No 1907/2006

mega TRAY NF

Revision date: 14.05.2018

Product code: 21371_EN

Page 3 of 7

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid dust formation. Do not breathe dust.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

6.4. Reference to other sections

Safe handling: see section 7 Personal protection equipment: see section 8 Disposal: see section 13

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

No special measures are necessary.

Advice on protection against fire and explosion

No special fire protection measures are necessary.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container tightly closed.

Advice on storage compatibility

No special measures are necessary.

7.3. Specific end use(s)

Production from individual trays

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Exposure limits (EH40)

| CAS No | Substance | ppm | mg/m³ | fibres/ml | Category | Origin |
|---------|--------------------|-----|-------|-----------|---------------|--------|
| 94-36-0 | Dibenzoyl peroxide | - | 5 | | TWA (8 h) | WEL |
| | | - | - | | STEL (15 min) | WEL |

8.2. Exposure controls

Protective and hygiene measures

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

Eye/face protection

Wear eye protection/face protection.

Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

Skin protection

Wear suitable protective clothing.



according to Regulation (EC) No 1907/2006

mega TRAY NF

Revision date: 14.05.2018

Product code: 21371_EN

Page 4 of 7

megadental GmbH

Respiratory protection

In case of inadequate ventilation wear respiratory protection.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

| Physical state: | Powder |
|--|--|
| Colour: | See packaging. |
| Odour: | characteristic |
| pH-Value: | not determined |
| Changes in the physical state | |
| Melting point: | ca. 110 °C |
| Initial boiling point and boiling range: | not determined |
| Flash point: | > 250 °C |
| Flammability | |
| Solid: | not determined |
| Gas: | not applicable |
| Lower explosion limits: | not determined |
| Upper explosion limits: | not determined |
| Ignition temperature: | > 400 °C |
| Auto-ignition temperature | |
| Solid: Gas: | not determined not applicable |
| Decomposition temperature: | not determined |
| | not determined |
| Oxidizing properties Not oxidising. | |
| Vapour pressure: | not determined |
| Density (at 20 °C): Bulk density (at 20 °C): | 1,16 g/cm³ 700 - 750 kg/m³ |
| Water solubility: | The study does not need to be conducted because the substance is known to be insoluble in water. |
| Solubility in other solvents Soluble in hydrocarbons. | |
| Partition coefficient: | not determined |
| Vapour density: | not determined |
| Evaporation rate: | not determined |
| 9.2. Other information | |
| Solid content: | 100 % |
| SECTION 10: Stability and reactivity | |

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

10.3. Possibility of hazardous reactions



according to Regulation (EC) No 1907/2006

mega TRAY NF

Revision date: 14.05.2018

Product code: 21371_EN

Page 5 of 7

No known hazardous reactions.

10.4. Conditions to avoid

none

10.5. Incompatible materials

No information available.

10.6. Hazardous decomposition products

Exothermal decomposition with formation of: Methyl methacrylate

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Toxicocinetics, metabolism and distribution

The product has not been tested.

Acute toxicity

Based on available data, the classification criteria are not met.

Irritation and corrosivity

Based on available data, the classification criteria are not met.

Sensitising effects

May cause an allergic skin reaction. (dibenzoyl peroxide; benzoyl peroxide)

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Additional information on tests

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

SECTION 12: Ecological information

12.1. Toxicity

The product is not: Ecotoxic.

12.2. Persistence and degradability

The product has not been tested.

12.3. Bioaccumulative potential

The product has not been tested.

12.4. Mobility in soil

The product has not been tested.

12.5. Results of PBT and vPvB assessment

The product has not been tested.

12.6. Other adverse effects

No information available.

Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

SECTION 13: Disposal considerations

13.1. Waste treatment methods



according to Regulation (EC) No 1907/2006

mega TRAY NF

Product code: 21371_EN

Page 6 of 7

Advice on disposal

Revision date: 14.05.2018

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information

Land transport (ADR/RID)

| <u>14.1. UN number:</u> | No dangerous good in sense of this transport regulation. |
|---|---|
| 14.2. UN proper shipping name: | No dangerous good in sense of this transport regulation. |
| <u>14.3. Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| 14.4. Packing group: | No dangerous good in sense of this transport regulation. |
| Inland waterways transport (ADN) | |
| <u>14.1. UN number:</u> | No dangerous good in sense of this transport regulation. |
| 14.2. UN proper shipping name: | No dangerous good in sense of this transport regulation. |
| <u>14.3. Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| 14.4. Packing group: | No dangerous good in sense of this transport regulation. |
| Marine transport (IMDG) | |
| <u>14.1. UN number:</u> | No dangerous good in sense of this transport regulation. |
| 14.2. UN proper shipping name: | No dangerous good in sense of this transport regulation. |
| 14.3. Transport hazard class(es): | No dangerous good in sense of this transport regulation. |
| 14.4. Packing group: | No dangerous good in sense of this transport regulation. |
| Air transport (ICAO-TI/IATA-DGR) | |
| <u>14.1. UN number:</u> | No dangerous good in sense of this transport regulation. |
| 14.2. UN proper shipping name: | No dangerous good in sense of this transport regulation. |
| <u>14.3. Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| 14.4. Packing group: | No dangerous good in sense of this transport regulation. |
| 14.5. Environmental hazards | |
| ENVIRONMENTALLY HAZARDOUS: | no |
| 14.6. Special precautions for user No information available. | |
| 14.7. Transport in bulk according to Anne not applicable | ex II of Marpol and the IBC Code |
| SECTION 15: Regulatory information | |
| 15.1. Safety, health and environmental re | gulations/legislation specific for the substance or mixture |

| EU regulatory information | |
|--|--|
| Information according to 2012/18/EU Not subject to 2012/18/EU (SEVESO III) (SEVESO III): | |
| Additional information | |
| To follow: 850/2004/EC, 79/117/EEC, 689/2008/EC | |
| National regulatory information | |
| Water contaminating class (D): 2 - clearly water contaminating | |
| Skin resorption/Sensitization: Causes allergic hypersensitivity reactions. | |



according to Regulation (EC) No 1907/2006

mega TRAY NF

Revision date: 14.05.2018

Product code: 21371_EN

Page 7 of 7

15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

SECTION 16: Other information

Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route

(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

Classification for mixtures and used evaluation method according to Regulation (EC) No. 1272/2008 [CLP]

| Classification | Classification procedure |
|--------------------|--------------------------|
| Skin Sens. 1; H317 | Calculation method |

Relevant H and EUH statements (number and full text)

Heating may cause a fire or explosion. May cause an allergic skin reaction. Causes serious eye irritation.

Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)